

Datasheet

Catalog Number: NBL-250779

Version: RN6.0




Product Name	NebuSelect™ Recombinant Human GDF9
Alias/Synonyms	GDF-9 Protein, zgc:123030 Protein, Gdf-9 Protein, growth differentiation factor 9 Protein, GDF9 Protein, gdf9 Protein, Gdf9 Protein
Appearance	Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.
Description	NebuSelect™ Recombinant Human GDF9, His-tag(Cat#NBL-250779) is expressed in Escherichia coli (E. coli) with His-tag at N-terminus.
Accession	O60383
AA Sequence	AlaGly320~Arg454
MW	The protein has a predicted MW of 16.8kDa.
Purity	> 95%
Endotoxin	< 0.1EU/ug(LAL Method)
Size	100ul; 500ul; 1ml
Storage	Short term 4°C; Long term -20°C. Avoid repeated freeze-thaw cycles.
Application	Positive Control; Immunogen; SDS-PAGE; WB.
Note	For research use only.

Figure 1 SDS-PAGE

Datasheet

Catalog Number: NBL-250779

Version: RN6.0



M: Marker; Lane 1: NBL-250779

4-20% SDS-PAGE Coomassie Staining (2ug)